INSTRUCTIONS FOR USE

Premium Teeth Resin

EN	Instructions for Use	NO	Instruksjoner for bruk
DE	Gebrauchsanweisung	RO	Instrucțiuni de utilizare
FR	Mode d'emploi	TR	Kullanım talimatları
IT	Istruzioni per l'uso	LV	Lietošanas instrukcija
ES	Instrucciones de uso	LT	Naudojimo instrukcijos
PT	Instruções de uso	ET	Kasutamise juhised
NL	Gebruiksaanwijzing	SK	Návod na použitie
DA	Brugsanvisning	SL	Navodila za uporabo
SV	Bruksanvisning	PL	Instrukcja użytkowania
FI	Käyttöohjeet	BG	Инструкции за употреба
CS	Návod k použití	KO	사용 지침
HU	Használati útmutató	JA	使用説明書
EL	Οδηγίες χρήσης	ZH	使用说明

Prepared: 10/29/2024 PRNT-0126 Rev 02

INSTRUCTIONS FOR USE

Premium Teeth Resin

Keep away from sunlight

Use-by date

Consult Instructions for Use

European Conformity

Batch Code

Catalog Number

Manufacturer

Temperature Limit

Prescription only

Medical Device

Environmental hazard

Warning: Read Safety Information

Unique device identifier

Importer

Indicates the authorized representative in Switzerland

Authorized representative in the European Community

Formlabs Ohio Inc.
27800 Lemoyne Rd,
Millbury, OH 43447
USA
+1 617 855 0762

EU Importer
Formlabs GmbH
Mühlenstraße 15, Floor 1
10243 Berlin, Germany
+49 3091 734306

Curaden AG
Amlehnstrasse 22
6010 Kriens
Switzerland

Advena Ltd.
Tower Business Centre
2nd Flr. Tower Street,
Swatar, BKR 4013 Malta

Before printing, refer to the Manufacturing Guide and dental applications guide for this material at support.formlabs.com for printing and post-processing recommendations and requirements to ensure the correct and safe usage of this material. For detailed safety and environmental information, the Safety Data Sheet is available at support.formlabs.com.

I. Introduction And Indications For Use

INTRODUCTION AND INTENDED PURPOSE

Premium Teeth Resin is a light-curable nano-ceramic filled resin designed for the fabrication of 3D printed dental appliances, such as denture teeth for complete and partial removable dentures, try-in dentures, provisional full-arch implant-supported restorations, and provisional restorations such as temporary crowns and bridges, inlays, onlays and veneers.

INDICATIONS FOR USE

Premium Teeth Resin when utilized to print 3D printed dental appliances, such as denture teeth for complete and partial removable dentures, try-in dentures, provisional full-arch implant-supported restorations, and provisional restorations such as temporary crowns and bridges, inlays, onlays and veneers is indicated to replace or restore missing tooth structures or missing teeth.

CONTRAINDICATIONS

Premium Teeth Resin should not be used in case of any known allergies or hypersensitivity history to photopolymer resin materials. In case of an allergic reaction, stop use immediately, and consult a physician.

INTENDED PATIENT POPULATION

All patients and populations in which a clinician has prescribed a treatment associated with the intended purpose.

INTENDED USERS

Dental and prosthodontic professionals who fabricate dental and prosthodontic appliances.

CLINICAL BENEFIT

The clinical benefit of Premium Teeth Resin is the ability to 3D print dental appliances, such as denture teeth for complete and partial removable dentures, try-in dentures, provisional full-arch implant-supported restorations, and provisional restorations such as temporary crowns and bridges, inlays, onlays and veneers.

II. Safety

RESIN ASSOCIATED WARNINGS

Premium Teeth Resin (uncured) contains polymerizable monomers.

- Please read the Premium Teeth Resin SDS at support.formlabs.com for hazards associated with liquid resin.

RESIN USE ASSOCIATED PRECAUTIONS

- Protective glasses and nitrile gloves should be worn while handling resin. For information about handling Premium Teeth Resin, please read the SDS at support.formlabs.com.

INFORMATION FOR PATIENT SAFETY

- The patient should inform the clinician of any known allergy or hypersensitivity history to photopolymer resins.
- In case of an allergic reaction, stop use immediately, and consult a medical physician.
- Do not use appliances if damage or cracks are present.
- If signs of extreme wear, deterioration, and/or perforations appear, the appliances should be replaced to avoid material failure and/or potential bodily injury.
- Device lifetime will vary from patient to patient, however, removable devices should be replaced at least every 5 years, and for provisional full-arch implant-supported restorations and provisional restorations including temporary crowns and bridges, inlays, onlays and veneers, devices should be replaced after a maximum period of 12 months.
- To maximize the device's lifetime, patients should be instructed to handle devices carefully and follow the cleaning and storage instructions.

III. Performance Characteristics

For detailed information on the mechanical properties of this material, technical data sheets can be found at dental.formlabs.com. Premium Teeth Resin has been evaluated in accordance with ISO 10993-1:2020, and ISO 7405:2018 and meets the requirements of all applicable biocompatibility testing.

IV. Specific Manufacturing Considerations

Refer to the Manufacturing Guide at support.formlabs.com for printing recommendations and requirements to ensure the correct and safe usage of this material.

- a. Hardware: Formlabs 3D Printers and accessories
- b. Software: Formlabs PreForm
- c. Part Orientation: Orient intaglio surface facing away from the build platform at a 0 - 30° tilted angle and generate supports automatically or manually.
- d. Part Geometry:
 - Single units (crowns, inlays, onlays, veneers):
 - Minimum wall thickness: 1.0 mm
 - Bridges (non full-arch):
 - [Anterior] Minimum wall thickness: 1.0 mm
 - [Anterior] Minimum cross sectional area of the connector: 12 mm²
 - [Posterior] Minimum wall thickness: 1.0 mm
 - [Posterior] Minimum cross sectional area of the connector: 16 mm²
 - Maximum construction length: 7 units
 - Pontic maximum: one molar width
 - Temporary Full-arch Implant-supported restorations (All-on-X):
 - Minimum wall thickness: 1.0 mm
- e. Recommended Post-Processing Equipment: Guidelines for Post-Processing Equipment can be found in the Manufacturing Guide.

REQUIREMENTS

For biocompatibility compliance, Premium Teeth Resin was validated with a dedicated resin tank, mixer, build platform, wash unit, and post-processing equipment that were not mixed with any other resins.

A. PRINTING

Printing guidelines and validated equipment can be found in the Manufacturing Guide at support.formlabs.com.

B. WASHING

Washing guidelines and validated equipment can be found in the Manufacturing Guide at support.formlabs.com.

C. POST-CURING

Post-Curing guidelines and validated equipment can be found in the Manufacturing Guide at support.formlabs.com.

D. SUPPORT REMOVAL

Support marks can cause abrasion if not removed and polished. Remove supports and raft using a cutting disk and handpiece, cutting plier, or other appropriate finishing tool. See Manufacturing Guide at support.formlabs.com for polishing recommendations.

E. CLEANING & DISINFECTION

1. Appliances may be cleaned using neutral soap and water or effervescent dental appliance cleaning tablets (used according to the manufacturer's directions).
2. Appliances may be disinfected by soaking in 70% concentration IPA for five minutes per CDC guidance (Guideline for Disinfection and Sterilization in Healthcare Facilities (2008)).
3. Inspect appliances for cracks after cleaning or disinfection. Discard if any damage or cracks are detected.

F. STORAGE

1. Printed appliances should be stored in closed, opaque or amber containers. Excess light exposure over time may affect the color and performance of printed appliances.
2. Store cartridges and covered resin tanks with printing resin in a dry place at 10°C - 25°C (50°F - 77°F). Do not exceed 25°C (77°F) when in storage.
3. Keep cartridges closed and away from ignition sources.
4. When transporting uncovered resin tanks, make sure that the printing resin is exposed to daylight for as short a period of time as possible.
5. Do not use resin past the expiration date listed on the cartridge label.

Notice: Expected results cannot be guaranteed if materials which have exceeded their minimum shelf life date are used or if storage instructions are not followed.

G. DISPOSAL

1. Cured resin is non-hazardous and may be disposed of as regular waste.
2. Follow facility protocols for waste that may be considered biohazardous.
3. Unused liquid resin & solvents should be disposed of in accordance with local regulations.
4. Cartridge & contaminated packaging should be disposed of in accordance with local regulations.
 - Contact a licensed professional waste disposal service to dispose of liquid resin.
 - Do not allow waste to enter storm or sewer drainage systems.
 - Avoid release into the environment.
 - Contaminated packaging: Dispose of cartridge or any unused product per local waste management guidelines.

H. ADVERSE EVENTS

To report any serious injury (temporary or permanent serious deterioration of health, a serious public health threat, or death), contact Formlabs Inc. (+1 617 855 0762) or your local Competent Authority.

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