VITA

VITA VIONIC® BASE RESIN FLEX

DE	Gebrauchsanweisung Kurzversion	4
EN	Instructions for use Short version see page 119 for US version	8
FR	Mode d'emploi Version courte	12
ES	Instrucciones de uso Versión abreviada	16
IT	Istruzioni d'uso Versione breve	20
BG	Инструкция за употреба Съкратена версия	24
CS	Návod k použití Stručná verze	28
DA	Brugsanvisning kort version	32
EL	Οδηγίες χρήσης Συνομειωμένη έκδοση	36
ET	Kasutusjuhend lühiversioon	40
FI	Käyttöohje Lyhyt versio	44
HR	Upute za uporabu kratka verzija	48
HU	Használati utasítás Rövid változat	52
JA	使用説明書 ショートバージョン	56
LT	Naudojimo instrukcija Trumpa versija	60

LV	Lietošanas instrukcija Īsā versija	64
NO	Bruksanvisning Kortversjon	68
NL	Gebruiksaanwijzing Korte versie	72
PL	Instrukcja używania Skrócona wersja	76
pt-BR	Instruções de uso Versão resumida	80
pt-PT	Instruções de utilização Versão resumida	84
RO	Instrucțiuni de utilizare Versiunea pe scurt	88
RU	Инструкция по применению Краткая версия	92
SK	Návod na použitie skrátená verzia	96
SL	Navodila za uporabo Kratka različica	100
SV	Bruksanvisning Kort version	104
TR	Kullanım kılavuzu Kısa versiyon	108
UA	Інструкція із застосування Коротка версія	112
ZH	使用說明書 簡短版	116
en-US	Instructions for use Short version	119

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VITA VIONIC BASE RESIN FLEX is a light-curing resin for the fabrication of flexible partial denture bases using a 3D printing process.

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Build process:

1. Prepare data:
 - CAD design.
 - Recommended **adhesive gap setting:**
 - VITA VIONIC BASE RESIN FLEX in combination with VITA VIONIC DENT RESIN or VITA VIONIC DENT DISC multiColor: **150 µm**
 - Generate an STL file.
 - Transfer the data to the CAM software and create support structures. Edit the support structures suggested by the software. Avoid support structures in the cavities.
 - Convert to a printer-compatible format.
2. Select process and material parameters in the printer software: Use the print template for VITA VIONIC BASE RESIN FLEX published by your printer manufacturer.

3. Transfer prepared data to the 3D printer.
4. Prepare 3D printing: Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller to avoid color deviations and misprints.
5. Fill the resin tank of the 3D printer according to the printer manufacturer's instructions.
6. Start the printing process.

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Final cleaning process:

15. Cleaning by the dental laboratory / practice is carried out in an ultrasonic bath. Do not use any thermal disinfection or sterilization methods. This could lead to deformation of the denture.
16. For daily cleaning of the denture by the patient, lukewarm water, a soft toothbrush, liquid soap or dish-washing detergent should be used. Rinse well afterwards. Denture cleaning tablets can be used according to the manufacturer's instructions.

Notes: Do not use toothpaste; the abrasive particles are coarse and scratch the denture.
 Do not use any acids or solvents.

General Specifications

- Follow the hardware manufacturer's instructions regarding parameter settings and printing and post-curing recommendations.
- A layer thickness of 100 µm is suitable for the printing process.
- The ideal processing temperature is 23 °C ± 2 °C / 70–77 °F.
- To avoid adverse effects on material quality, do not expose the liquid material to light or UV irradiation under any circumstances.
- Deviations from the described manufacturing processes or storage conditions can lead to deviating mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.

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The devices should be set up, validated and maintained in accordance with their labeling and instructions for use.

An up-to-date list of **validated 3D printers, washing units and post-curing units** as well as suitable printer templates and post-curing parameters can be found at www.vita-zahnfabrik.com/VionicBaseResinFlex.

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Safety Information

- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- For open handling, use devices with local exhaust ventilation.
- Wear suitable protective equipment when working (when handling the uncured material: Protective gloves and safety goggles, when post-processing the cured material, a face mask should also be worn).

- Avoid contact with the skin. In case of contact with skin, wash immediately with plenty of soap and water.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Biocompatibility is guaranteed only after complete polymerization.
- MMA-free.
- Refer to the relevant safety data sheet for hazard and safety information.
- Any serious incident that has occurred in connection with this product must be reported immediately to info@vita-zahnfabrik.com and to the competent public authority of the Member State in which the user and/or patient resides.

Disclaimers

Please note: Our products must be used in accordance with the instructions for use. We accept no liability for damage resulting from improper handling or processing. The user shall also test the product for its suitability for the intended application before using it. Liability on our part is excluded if the product is processed in a non-contractual or non-permissible combination with materials and/or devices from other manufacturers and this results in damage. Publication of these Instructions for Use: 2025-12.

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Instructions for Use I Short Version

Please read the Instructions for Use before using the product.

VITA VIONIC BASE RESIN FLEX is a light-curing resin for the fabrication of flexible partial removable denture bases using a 3D printing process. It can be used in combination with a stereolithographic 3D printer using a 385 nm light source. The product processing requires a computer aided and manufacturing system (CAD/CAM) that includes the following: scanner, design software, additive printer and post-cure unit. Those devices and software are not part of the product.

CAUTION STATEMENT: Caution: Federal Law restricts this device to sale by or on the order of a dentist.

Intended Use

Resin for dental 3D printing.

Patient Target Group

Persons being treated in the context of a dental procedure.

Intended Users

Exclusively specialist users: Dentists and dental technicians.

Indication

For the fabrication of flexible denture bases.

Contraindication

Contains (meth)acrylates and phosphine oxides.

Ingredients of VITA VIONIC BASE RESIN FLEX may cause allergic reactions in accordingly predisposed persons. In such a case, do not use the product any further. Apply the product intraorally only when it is fully polymerized.

Product Application/Operation

Please observe the following workflow during processing:

Build process:

1. Prepare data:

- CAD design.

Recommended **adhesive gap setting:**

- VITA VIONIC BASE RESIN FLEX in combination with VITA VIONIC DENT RESIN or VITA VIONIC DENT DISC multiColor: **150 µm**

- Generate an STL file.

- Transfer the data to the CAM software and create support structures. Edit the support structures suggested by the software. Avoid support structures in the cavities.
 - Convert to a printer-compatible format.
2. Select process and material parameters in the printer software: Use the print template for VITA VIONIC BASE RESIN FLEX published by your printer manufacturer.
 3. Transfer prepared data to the 3D printer.
 4. Prepare 3D printing: Before use, the material in the bottle should be shaken intensively or homogenized with a bottle roller to avoid color deviations and misprints.
 5. Fill the resin tank of the 3D printer according to the printer manufacturer's instructions.
 6. Start the printing process.

Post-processing:

7. Allow the printed objects to drip dry in the printer for approx. 10 minutes. Ideally, the post-processing should take place immediately after the build process.
8. Take the platform out of the printer and remove the printed objects from the platform. If there are several objects, please separate them.
9. Remove the outer support structures.
10. Remove excess resin with compressed air.

Cleaning:

11. Clean the printed objects twice for 3 minutes each in an ultrasonic bath, using fresh isopropanol ($\geq 98\%$) each time, and follow the device manufacturer's instructions.
12. Dry the cleaned print objects until there are no isopropanol residues left. Check openings, cavities and gap areas for residues, and blow them off with compressed air for faster drying. After cleaning, the surface must no longer be sticky and shiny.

Post-exposure:

13. Carry out post-curing in a suitable light-curing device in accordance with the manufacturer's instructions. The properties and color appearance of the final product depend, among other things, on post-processing. Correct post-exposure is essential for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured.

Surface processing:

14. The denture base is now ready for further processing or finishing. Mechanically process the surface (smooth with a milling cutter) and polish. Pre-polishing is done with brushes and pre-polishing paste/pumice, high-gloss polishing with buffing and high-gloss polishing agents for resins.

Final cleaning process:

15. Cleaning by the dental laboratory / practice is carried out in an ultrasonic bath. Do not use any thermal disinfection or sterilization methods. This could lead to deformation of the denture.
16. For daily cleaning of the denture by the patient, lukewarm water, a soft toothbrush, liquid soap or dish-washing detergent should be used. Rinse well afterwards. Denture cleaning tablets can be used according to the

manufacturer's instructions.

Notes: Do not use toothpaste; the abrasive particles are coarse and scratch the denture.
Do not use any acids or solvents.

General Specifications

- Follow the hardware manufacturer's instructions regarding parameter settings and printing and post-curing recommendations.
- A layer thickness of 100 µm is suitable for the printing process.
- The ideal processing temperature is 23 °C ± 2 °C / 70–77 °F.
- To avoid adverse effects on material quality, do not expose the liquid material to light or UV irradiation under any circumstances.
- Deviations from the described manufacturing processes or storage conditions can lead to deviating mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.
- It is not allowed to pour 3D printing resin from the reservoir of the printer back into the bottle.

Recommended Product Combinations

For recommended product combinations with the VITA system components as well as corresponding processing recommendations and restrictions regarding the combinability of the product, please refer to the Instructions for Use | Full Version.

Use of the product is permissible only with compatible devices.

The devices should be set up, validated and maintained in accordance with their labeling and instructions for use. An up-to-date list of **validated 3D printers**, **washing units** and **post-curing units** as well as suitable printer templates and post-curing parameters can be found at www.vita-zahnfabrik.com/VionicBaseResinFlex.

Storage/Disposal

Storage at room temperature 15 °C (59 °F) to 28 °C (82 °F). Protect from direct sunlight. Even slight exposure to light can trigger polymerization. To protect from contamination, cover the material in the tray with the lid or a glass plate. Products marked with a hazardous-substance pictogram must be disposed of as hazardous waste. Recyclable waste (such as attachments, paper, resins and plastics) must be disposed of via appropriate recycling systems. If necessary, contaminated product residues must be pretreated and disposed of separately in accordance with regional/national and international regulations and provisions.

Adverse Effects

Product may cause allergic reactions.

Safety Information

- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- For open handling, use devices with local exhaust ventilation.
- Wear suitable protective equipment when working (when handling the uncured material: Protective gloves and safety goggles, when post-processing the cured material, a face mask should also be worn).
- Avoid contact with the skin. In case of contact with skin, wash immediately with plenty of soap and water.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Biocompatibility is guaranteed only after complete polymerization.
- MMA-free.
- Refer to the relevant safety data sheet for hazard and safety information.
- Any serious incident that has occurred in connection with this product must be reported immediately to info@vita-zahnfabrik.com and to the competent public authority of the Member State in which the user and/or patient resides.

MRI Safety Information: VITA VIONIC BASE RESIN FLEX is MR Safe as it is composed of materials that are electrically non conductive, nonmetallic and nonmagnetic.

Technical Data

Chemical composition	Wt%
Urethane dimethacrylate oligomer	40-60 %
Methacrylate monomers	15-60 %
Fillers, pigments, stabilizers	< 5 %
Initiators	0.2-5 %
MMA and phthalat free	0%

Physical properties	Value
Flexural strength DIN EN ISO 20795-1*	≥ 25 MPa
Flexural modulus DIN EN ISO 20795-1*	≥ 900 MPa
Fracture work DIN EN ISO 20795-1*	> 7350 J/m ²
Water sorption DIN EN ISO 20795-1*	< 32 µg/mm ³
Water solubility DIN EN ISO 20795-1*	<3.5 µg/mm ³
Viscosity	1000-1700 mPas

*Applies to cured objects | according to norm at room temperature

Disclaimers

Please note: Our products must be used in accordance with the instructions for use. We accept no liability for damage resulting from improper handling or processing. The user shall also test the product for its suitability for

the intended application before using it. Liability on our part is excluded if the product is processed in a non-contractual or non-permissible combination with materials and/or devices from other manufacturers and this results in damage. Publication of these Instructions for Use: 2025-12.

With the publication of these instructions for use, all previous editions lose their validity. You can find the latest version at www.vita-zahnfabrik.com.

Symbol Explanation

Graphic	Title	Description	Reference
	Manufacturer	Indicates the medical device manufacturer.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Distributor	Indicates the medical device manufacturer.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Medical device	Indicates that the product is a medical device.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Prescription only	Requires a prescription in the United States of America. Caution: In the United States of America, federal law restricts this device to sale or use by, or on the order of, a physician.	USA Code of Federal Regulations 21CFR Part 801 § 801.109 (b)(1)

	Magnetic resonance	The freeprint crown is MR Safe.	FDA Guidance document: Testing and Labeling Medical Devices for Safety in the Magnetic Resonance (MR) Environment dated May 20, 2021.
	Caution	The uncured material is considered as hazardous substance.	29 CFR 1910.1200 Occupational Safety and Health Standards
	This way up	This is the correct upright position of the distribution packages for transport and/or storage.	ISO 780:2015(E) Packaging-Distribution packaging- Graphical symbols for handling and storage of packages

VITA Zahnfabrik has been certified, and the product **VITA VIONIC®** BASE RESIN FLEX bears the mark:

CE 0483

Made in Germany

VITA Zahnfabrik H. Rauter GmbH & Co. KG

Spitalgasse 3
79713 Bad Säckingen
Germany

Phone: +49 7761 562-0
Fax: +49 7761 562-299
E-Mail: info@vita-zahnfabrik.com

detax GmbH

Carl-Zeiss-Str. 4
76275 Ettlingen
Germany

Phone: +49 7243 510-0
Fax: +49 7243 510-100
E-Mail: post@detax.com

VITA